



AN INFORMATIONAL RESOURCE FOR PROVIDERS OF AND SUPPLIERS TO EMS

Treating Congestive Heart Failure with CPAP

Congestive Heart Failure (CHF) is defined as, “excessive blood or fluid in the lungs or body tissues caused by the failure of ventricles to pump blood effectively”.

The condition is termed congestive because the fluids congest, or clog the organs. It is termed heart failure because the congestion both results from and also aggravates failure of the heart to function properly.

Congestive Heart failure occurs when the left ventricle cannot pump out the amount of blood entering the ventricle, or when the ventricle is damaged and cannot effectively pump enough blood to meet the body’s requirements. It may also occur due to a build up of excess fluid in the body due to kidneys being damaged or not functioning properly due to disease.

Blood begins to congest in the lungs (pulmonary edema). The work of breathing increases as the airways are obstructed by the fluids, impeding the flow of air into the lungs. The alveoli are unable to exchange gases effectively creating severe shortness of breath in patients. As the condition worsens, this congestion will eventually cause the right ventricle to fail. When this occurs, valuable blood supply (oxygen and nutrients) to the body’s cells are seriously disrupted, and the waste products of metabolism are no longer eliminated effectively causing these toxins to accumulate and ultimately causing cell death.

CPAP (Continuous Positive Airway Pressure)

EMS Providers play a large role in the emergency management of the patients suffering from Congestive Heart Failure.

CPAP is a form of Noninvasive Positive Pressure Ventilation (NPPV) which is becoming increasingly popular in the field management of the patient suffering from CHF. CPAP can save precious seconds when managing a patient with Congestive Heart Failure, avoiding intubation and its associated field and long-term concerns and side effects.

CPAP improves the ability of the alveoli to diffuse oxygen to the red blood cells, by using pressure to drive gas into the alveoli and open up unused or collapsed alveoli. CPAP increases the resistance of gas flow during exhalation providing resistance to the exiting airflow of gas from the lungs.

CPAP can lessen the shortness of breath experienced by the patient by improving the performance of the heart. Ventilation improves and airway secretions are removed improving oxygenation and CO₂ removal.

CPAP also helps by increasing intrathoracic pressure, causing an increase in cardiac output. It helps to reduce the need for CHF patients to be intubated and placed on a ventilator.

The patients must be able to breathe spontaneously to benefit from this patient management protocol. The patient breathes through a pressurized circuit against a set threshold resistor that keeps the airway pressure at a preset level, normally monitored on an airway pressure gauge. The inspiratory and expiratory pressures remain the same.

If the patient moves into respiratory arrest during CPAP treatment, artificial ventilation will be required via a bag-valve-mask or automatic ventilator. Devices like the CAREvent® ALS+CPAP offer both CPAP and ventilation in one space saving design to treat patients moving from difficulty to not breathing.

Management of a patient suffering from Congestive Heart Failure in the field with CPAP, increases the chance that the patient can avoid an ICU stay on a ventilator and reducing the associated costs

Street Sense is published by O-Two Medical Technologies Inc., innovators in resuscitation since 1971 (mouth-to-mouth barriers, BVMs to a full line of automatic ventilators/resuscitators).

Contact O-Two at:

7575 Kimbel St.
Mississauga ON Canada L5S 1C8
Email: mboadway@otwo.com
Phone: 1-800-387-3405
Fax: 905-677-2035
Web : www.otwo.com