

Burn Relief®

Burn Management Products



Burn Relief® Large Burn Kit
Part #: 06BK1200



Burn Relief® Fire Blanket
Part #: 06BK1232



MEDICAL TECHNOLOGIES INC

"INNOVATION IN RESUSCITATION"

Burn Relief®

Burn Relief® helps to reduce trauma to your burn patients while simultaneously providing an effective “One Step” approach to the management of burns. Simply apply the **Burn Relief®** product to the burned area and the sterile, self-cooling, non-adhering **Burn Relief®** quickly helps to reduce trauma and burn progression for the burn patient while relieving pain.

The 96% water, gel-saturated, non-woven polyester composition of **Burn Relief®** draws heat away from the burn area. This provides the cooling effect which can last up to 8 hours!

As a solid, gelled water product, **Burn Relief®** is not wet or messy, has no pungent odour, does not absorb wound exudate, leaves no residue on the skin and contains no anti-microbial agents that may require flushing which could delay follow-up Emergency Room treatment.

Easily applied, **Burn Relief®** comes in a variety of sizes and configurations to meet the requirements of emergency burn management.

The large, 50” x 72” **Burn Relief® Fire Blanket** provides full body coverage for large area burns and also acts as an escape blanket. The non-woven polyester is fire resistant and the gel saturation allows the Fire Blanket to withstand temperatures up to 2000°C (3568°F) for several minutes. This makes the Fire Blanket a valuable, lifesaving rescue tool.

Burn Relief® FEATURES

- **Sterile:** Sterilized by radiation
- **Not Wet or Messy:** Solid gel does not drip or run
- **Self-Cooling:** Cools by evaporation and convection. Gel is 96% water and draws heat away from the burn area
- **Non-irritating:** Assures lack of irritation and further contamination
- **Non-Adhering:** The moist surface of the dressings prevents adherence to burned body tissue
- **Non-absorbent:** The gel is non-soluble in body fluids
- **Conforming:** Conforms to body areas for total contact without hindering patient movement
- **Protective:** Protects against further exposure of the wound to contamination

NOTE:

- Do not allow **Burn Relief®** Dressings to freeze. Store above 0°C (32°F) (**Burn Relief®** products are still effective after thawing).
- Sterility guaranteed providing packaging is not damaged or opened.
- Chemical burns must be thoroughly flushed before **Burn Relief®** is applied.
- Single use only.

Ordering Information:

06BK1200	Large Burn Kit - 1 x Body Towel 24" x 36", 2 x Face Mask 14" x 18", 2 x Leg Roll, 8" x 36", 2 x Arm Roll 3" x 36", 2 x Finger Roll 1" x 36", 2 x Self Conforming Gauze, 1 x Scissors, 1 x Carry Case
06BK1201	Mini Burn Kit - 1 x 8" x 36", 2 x 3" x 36", 2 x 1" x 36", 1 x Scissors, Carrying Case
06BK1203	Burn Pad 7.5 cm X 7.5 cm (3" X 3")
06BK1210	Burn Pad 10 cm X 10 cm (4" X 4")
06BK1220	Burn Roll 2.5 cm X 91.4 cm (1" X 36")
06BK1221	Burn Roll 7.5 cm X 91.4 cm (3" X 36")
06BK1222	Burn Roll 20.3 cm X 91.4 cm (8" X 36")
06BK1230	Face Mask 35.6 cm X 45.7cm (14" X 18")
06BK1231	Body Towel 61 cm X 91.4 cm (24" X 36")
06BK1232	Fire Blanket 127cm X 183 cm (50" X 72")