

CPR Personal Protective Devices

***“Bacterial & Viral Protection
During Training and CPR”***



Bio-Barrier® Faceshields



Revive Aid® Faceshields



Rescue Breather® CPR Pocket Ventilator



MEDICAL TECHNOLOGIES INC

“INNOVATION IN RESUSCITATION”

CPR Personal Protective Devices

The CPR Bio-Barrier® Faceshields

The 2-way, biological filter of the **CPR Bio-Barrier Faceshield** will prevent the passage of over 99% of the viruses and bacteria that cause diseases such as Tuberculosis, Hepatitis and HIV*! It also acts as a physical barrier to moisture, helping to prevent the transfer of body fluids.

During CPR training or any breathing/cardiac arrest emergencies the large, pliable, yet extremely strong faceshield, with graphic illustrations, assists the student / rescuer in correctly positioning the **CPR Bio-Barrier Faceshield**.

Packaged in a compact and durable "easy open" polyethylene pouch, the **CPR Bio-Barrier Faceshield** fits easily into wallets, purses and First Aid kits ready for any breathing emergency. It is also available in a convenient, attractive keychain pouch.

For training purposes the **CPR Bio-Barrier Faceshields** are packaged on a roll allowing each student to tear off their own personal protective device.

Features:

- Bacterial/Viral Filtration
- Clear Faceshield
- Faceshield is pliable yet strong
- Positioning illustrations

***Independent laboratory tested**

Revive Aid® CPR Protective Barrier

The single use, **Revive Aid® CPR Protective Barrier** is ready when you need it! Designed as a "First Response" device for use during CPR, the **Revive Aid®** provides a physical barrier between the rescuer and the patient.

This single use, cost-effective quality device provides rescuer protection by the use of a one-way flapper valve and filter to reduce the risk of cross contamination.

Compact and lightweight the **Revive Aid®** can be easily stored in First Aid kits, purses, sports bags and glove compartments. Packaged complete with instructions for use in a clear polyethylene Ziploc® bag, the **Revive Aid®** is also configured in a keychain pouch or key chain/belt pouch complete with a pair of gloves which can be worn by Policemen, Firefighters, Paramedics, First Responders, other Professional Healthcare Providers or anyone trained in CPR, ensuring that you will always carry the level of protection you deserve when helping to save a life.

Features:

- One size fits all patients.
- Ergonomic patient mouthpiece.
- Round rescuer mouth piece
- One-way valve with filter.
- Clear faceshield gives excellent coverage yet allows the rescuer to see the patient's face.

Rescue Breather® CPR Pocket Ventilator

The **Rescue Breather® CPR Pocket Ventilator** provides the rescuer with a safe and effective means to perform life saving expired air resuscitation while offering maximum protection from the potential risks of disease transfer for both the patient and the rescuer.

The **Rescue Breather® CPR Pocket Ventilator** comes standard with a low resistance disposable one way valve, with 3M Filtrete® filter, oxygen inlet, head strap and carrying case.

The filter ensures optimum bacterial protection, yet minimum airway resistance. The seamless air cushion mask eliminates seam failures and air loss during storage.

The compact size of the **Rescue Breather® CPR Pocket Ventilator** makes it able to fit into your pocket, purse, brief case, sports bag or glove compartment of your car, keeping it ready for any breathing emergency.

Features:

- Transparent Universal Cushion Mask for Adults to Infants
- Replaceable, clear, one-way valve with 3M Filtrete® filter
- Oxygen Inlet with protective cap
- Head Strap
- Compact, "Clamshell" carrying case with hold tabs

Ordering Information:

01RB5010 **CPR Bio-Barrier® Faceshield** in polyethylene pouch (200/case)
01RB5030-cs **CPR Bio-Barrier® Faceshield** in keyring pouch (50/case)
01RB5050 **CPR Bio-Barrier® Training Shield** in dispenser roll of 50 (4 rolls per case)

Ordering Information:

01RB1000-cs **Revive Aid® CPR Device** in clear polyethylene bag (cs/50)
01RB2010Red **Revive Aid® CPR Device** in Red nylon keychain pouch (cs/50)
01RB2010Black **Revive Aid® CPR Device** in Black nylon keychain pouch (cs/50)
01RB2020-cs **Revive Aid® CPR Device** in nylon keychain/belt pouch c/w gloves (cs/50)
(Available in Red, Black and Blue)

Ordering Information:

01PV1500W-cs **Rescue Breather® CPR Pocket Ventilator** c/w Headstrap, Oxygen Inlet & One-Way Valve in a Clam Shell Hard Case with 3M Filtrete® Filter - Case / 12
01PV1505W-cs **Replacement One-Way Valves** - Pkg./ 10

NOTE: All devices shown are for emergency use only by persons trained in CPR and in the use of these types of devices. The use of these devices does not guarantee 100% protection against disease transmission