

Easy Grip[®]

Bag-Valve-Mask Resuscitators

“Unbelievable Quality at an Unbelievable Price”!



Easy Grip[®] Single Patient Use



“INNOVATION IN RESUSCITATION”

Easy Grip®

Bag-Valve-Mask Resuscitators

The **Easy Grip**® Bag-Valve-Mask Resuscitators are manual resuscitation devices designed for use in various clinical settings to provide respiratory support to the respiratory distressed or non-breathing patient.

The Child and Infant Resuscitators come complete with a pressure Relief Valve set at 40 cm H₂O, that can be overridden by the operator or locked into position.

When used as an oxygen supplemental device, these devices can provide up to 100% oxygen on every squeeze and release of the bag.

Available in three versions:

**Disposable
Autoclavable Silicone
and Cold Chemical Sterilizable**

and three sizes:

**Adult
Child
Infant**

The **Easy Grip**® range offers the user a complete choice of the size and type of BVM to suit their patient and institutional needs.

NOTE: An optional pressure relief system is available for the Adult BVM (relief setting 60 cmH₂O)

CAUTION: Federal (U.S.A.) law restricts this device to sale by, or on the order of a physician.

WARNINGS:

Easy Grip® Bag-Valve-Mask Resuscitators should only be used by trained personnel familiar with Cardio-Pulmonary Resuscitation. Use only as directed. Improper usage or unauthorized modification of this product may result in user or patient injury.

Easy Grip® Disposable Bag-Valve-Mask Resuscitators are designed for single patient use only. They must not be cleaned and re-used.

SPECIFICATIONS

Storage Temperature:	-40°C to 60°C (-40°F to 140°F)	Deadspace:	7.0 ml
Operating Temperature:	-18°C to 50°C (0°F to 122°F)	Reservoir Volume:	2700 ml
Inspiratory Resistance:	3.3cm H ₂ O	Expiratory Resistance:	2.2cm H ₂ O

Ordering information:

O1BM3500-Cs	Disposable Easy Grip ® (Adult) c/w facemask, reservoir system and oxygen tubing - case/12
O1BM3510-Cs	Disposable Easy Grip ® (Child) c/w facemask, reservoir system and oxygen tubing - case/12
O1BM3520-Cs	Disposable Easy Grip ® (infant) c/w facemask, reservoir system and oxygen tubing - case/12
O1BM3600	Silicone Easy Grip ® (Adult) c/w facemask and reservoir system - each
O1BM3610	Silicone Easy Grip ® (Child) c/w facemask and reservoir system - each
O1BM3620	Silicone Easy Grip ® (infant) c/w facemask and reservoir system - each
O1BM3700	Cold Chemical Easy Grip ® (Adult) c/w facemask, reservoir system - each
O1BM3710	Cold Chemical Easy Grip ® (Child) c/w facemask, reservoir system - each
O1BM3720	Cold Chemical Easy Grip ® (infant) c/w facemask, reservoir system - each



O-TWO MEDICAL TECHNOLOGIES INC.
"INNOVATION IN RESUSCITATION"

7575 Kimbel Street, Mississauga, Ontario, Canada L5S 1C8
Toll Free: (800) 387-3405 Tel: (905) 677-9410 Fax: (905) 677-2035
E-mail: resuscitation@otwo.com
Website: www.otwo.com

