

Resuscitation & Anesthesia Facemasks



***“A Comprehensive Range
of Facemasks at the
Right Price”***

Resuscitation & Anesthesia Facemasks

With a wide choice in style, material and size, O-Two Medical Technologies Inc. provides a comprehensive range of resuscitation and anesthesia facemasks to meet any respiratory requirement.



From the truly Universal, “one-size-fits-all!” Facemask (available in both PVC and silicone) to the disposable, inflated cuff, range for adults to neonates, the high quality and superior design of these masks makes them the perfect choice for both hospital and pre-hospital use.

The silicone range includes “Donut” style neonate/infant masks as well as the flanged design of the O-Two Medical Technologies’ “Everseal Edge” range. These are ideal anesthesia masks as they will withstand autoclaving.

The three sizes of clear PVC, Everseal masks, are ideal for pre-hospital use where disinfection after use is normally undertaken by cold chemical sterilant such as an activated glutaraldehyde or a 10% bleach solution.

Where a single-use mask is required, the Disposable Inflated Cuff range meets the users needs with a range of three products with a soft sealing cushion and clear plastic body.



All masks have standard 22mm or 15mm (infant/neonate) inlet connectors and they come conveniently packaged in cases of 12 units.

Ordering Information:

- 02FM4999-cs **O-Two Medical Universal Mask** (PVC) - Case / 12
- 02FM5000-cs **O-Two Medical Medium Mask** (#2) with Anatomical Chin Cup (PVC) - Case / 12
- 02FM5001-cs **O-Two Medical Large Mask** (#3) with Anatomical Chin Cup (PVC) - Case / 12

- 02FM5097-cs **O-Two Medical Universal Mask** (Silicone) - Case / 12
- 02FM5098-cs **O-Two Medical Medium Mask** (#2) with Anatomical Chin Cup (Silicone) - Case / 12
- 02FM5099-cs **O-Two Medical Large Mask** (#3) with Anatomical Chin Cup (Silicone) - Case / 12
- 02FM5100-cs **O-Two Medical Toddler/Child Mask** (Silicone) - Case / 12
- 02FM5101-cs **O-Two Medical Infant Mask** (Silicone) - Case / 12
- 02FM5102-cs **O-Two Medical Neonate Mask** (Silicone) - Case / 12

- 02FM5205-cs **O-Two Medical Disposable Cuffed Neonate Mask** (#1) - Case / 12
- 02FM5207-cs **O-Two Medical Disposable Cuffed Child/Small Adult Mask** (#3) - Case / 12
- 02FM5210-cs **O-Two Medical Disposable Cuffed Adult Mask** (#4) - Case / 12



O-TWO MEDICAL TECHNOLOGIES INC.
“INNOVATION IN RESUSCITATION”

7575 Kimbel Street, Mississauga, Ontario, Canada L5S 1C8
Toll Free: (800) 387-3405 **Tel:** (905) 677-9410 **Fax:** (905) 677-2035
E-mail: resuscitation@otwo.com **Website:** www.otwo.com

