

OXYGEN EQUIPMENT CARRYING CASES AND AIRWAY MANAGEMENT KITS



Emergency Responder Moulded Case



Semi-Rigid Airway Management Case



Softpack Clamshell Case



Automatic Transport Ventilator Case

OXYGEN EQUIPMENT CASES AND AIRWAY MANAGEMENT KITS

Emergency Responder Moulded Plastic Carrying Case

A rugged performer that allows full access to your oxygen therapy and resuscitation equipment. Designed to accommodate a "D" or "Jumbo D" oxygen cylinder and Statvac II Oxygen-Powered Aspirator, the Emergency Responder is ready for any breathing emergency.



The Emergency Responder comes standard with a moulded transport hole in the end of the case to allow an oxygen mask or oxygen powered resuscitator to be used during transport with the lid of the case in a closed position.

Clamshell Carrying Case

The safe way to transport your oxygen therapy kit and supplies. The clamshell design allows the kit to be completely opened its entire length for easy access to all your emergency equipment. Additional storage for your respiratory accessories is easily accommodated by the large elastic top pouches along the inner wall. Large loop handles and adjustable shoulder strap allows hands free operation during transport.



Airway Management Kit

Semi-rigid oxygen therapy and resuscitation system that organizes your oxygen cylinder and regulator, resuscitator, oxygen powered aspirator and emergency airway management equipment. Large Top and side-opening main compartment with large accessory pouches that allows the lid to be fully opened for quick and easy access. The large Central compartment with adjustable oxygen cylinder straps to accommodate a range of cylinder sizes and styles. The end pocket has been specially designed to hold and operate an oxygen powered suction unit with accessories. Hidden adjustable backpack straps allow hands free operation during transport.



Automatic Transport Ventilator Case

The best way for you to store, protect and carry your **CAREvent**® series of Transport Ventilators while transporting a critically ill patient. Hands-free transport is easily accomplished by hanging the case over the rails of a stretcher or hospital bed. It features an easy open large zippered top and front lid. The lid also contains a removable accessory pouch that can accommodate your ventilator respiratory accessories.



Ordering information:

- 01CC8005** *Emergency Responder Moulded Carry Case* in International Orange for "D" or Jumbo "D" Cylinders (for oxygen therapy or resuscitation equipment)
- 01CC8016** *Softpack Clamshell Carrying Case* with sling Style Shoulder Strap for "D" or Jumbo "D" Cylinders (for oxygen therapy or resuscitation equipment)
- 01CC8025** *Airway Management Semi-Rigid Carrying Case* with "Hidden Backpack Shoulder Harness" and Sling Style Shoulder Strap (Designed to hold **CAREvent**® Automatic Transport Resuscitators and "D" or Jumbo "D" Cylinders)
- 01CC8030** *CAREvent*® *Automatic Transport Ventilator carrying case* with Sling Style Shoulder Strap and mounting hooks for stretcher attachment (specifically designed for **CAREvent**® control module units).



O-TWO MEDICAL TECHNOLOGIES INC.
"INNOVATION IN RESUSCITATION"

7575 Kimbel Street, Mississauga, Ontario, Canada L5S 1C8
Toll Free: (800) 387-3405 **Tel:** (905) 677-9410 **Fax:** (905) 677-2035
E-mail: resuscitation@otwo.com **Website:** www.otwo.com

