

Oxygen Therapy Kits



01OT3000 & 01OT3200 Inhalator Kits

**Oxygen Therapy
On The Move!**

Oxygen Therapy Kits

O-Two Medical Technologies Oxygen Therapy Kits are ready for any breathing emergency when seconds count. They can be ordered either in a **Emergency Responder Moulded Plastic Carrying Case** constructed in high visibility International Orange or our Royal Blue **Softpack Clamshell Carrying Case** with Sling Style Shoulder Strap for “hands-free” operation.

The O-Two Medical Technologies **SQ^T Brass Adjustable Oxygen Therapy Regulator** features high quality manufacturing and performance. The “Aluminum Free” brass construction improves product safety by reducing the risk of exacerbating regulator fires should they occur.

With its offset gauge and small diameter body, the **SQ^T** is one of the slimmest regulators of its type.

Ideal for all applications where adjustable oxygen therapy from 0 - 25 lpm flowrates are required. It is unbeatable for quality, size and performance!

The **Emergency Responder Moulded Plastic Carrying Case** has a built in carrying handle that allows the kit to be carried in a horizontal position and remains in a stable position during use. A moulded transport hole in the end of the case allows the oxygen therapy mask to be used with the case lid closed during patient transport. The case design makes the kit lightweight and buoyant should the kit fall into the water during an emergency. It will stay afloat for a short period of time allowing for immediate retrieval.

The **Softpack Clamshell Carrying Case** is constructed with extra thick closed cell foam, elastic airway holders, double reinforced abrasion resistant nylon bottom protecting your equipment from the rigors of Emergency use.

It features a large two way, self-repairing zipper that allow the case to be opened fully during any breathing emergency. Large, oversized, interior pockets allow the easy access storage of oxygen therapy disposables and accessories. The outside features reflective material for visibility and double carrying handles for security and an adjustable shoulder strap sling that frees the Rescuer hands for other tasks.

SPECIFICATIONS

SQ^T Oxygen Therapy Regulator

Therapy Flow Settings SQ ^T :	0, 0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25 (+/- 10%)
Gauge Accuracy:	2% of full scale +8% of reading
Maximum Rated Input Pressure:	2216 psi
Operating Temperature Range:	-0°F to +140°F
Storage Temperature Range:	-20°F to +140°F
Size SQ ^T (approx):	L 6" x H 1.75" x W 1.25"
Weight SQ ^T :	24.0 oz
Outlet Style:	Barbed nipple

Carrying Cases

Moulded Carry Case:

Size (approx):	L 23" x H 6.25" x W 10.25"(approx)
Color:	International Orange

Clamshell Softpack:

Size (approx):	L 19" x Dia. 8" (approx)
Color:	Royal Blue

Ordering Information:

- 01OT3000 O-Two Medical Inhalator Kit** c/w Adjustable 0 - 25 LPM Therapy Flow Regulator, Cylinder Handwheel, Medium Concentration Facemask with Tubing in a Moulded Carrying Case (Cylinder not Included)
- 01OT3200 O-Two Medical Inhalator Kit** c/w Adjustable 0 - 25 LPM Therapy Flow Regulator, Cylinder Handwheel, Medium Concentration Facemask with Tubing in a Softpack Clamshell Carrying Case (Cylinder not Included)



O-TWO MEDICAL TECHNOLOGIES INC.
“INNOVATION IN RESUSCITATION”

7575 Kimbel Street, Mississauga, Ontario, Canada L5S 1C8
Toll Free: (800) 387-3405 **Tel:** (905) 677-9410 **Fax:** (905) 677-2035
E-mail: resuscitation@otwo.com **Website:** www.otwo.com

